



Version 22.0 NEW FEATURES	Platform	ChemDraw Prime	ChemDraw Professional	ChemOffice	ChemOffice+ Cloud Standard*
Atom/Bond Grab Hotkey	Win/Mac	•	•	•	•
3D Rotation Shortcut	Win/Mac	•		•	
Open CIF Files	Win/Mac	•	•	•	•
Toggle CIF Files Hydrogens	Win/Mac	•	•	•	•
Hydrogen Bond Tool	Win/Mac				•
Millipore Sigma Synthia Add-in (requires Synthia subscription)	Win/Mac		•	•	•
HELM Cartoon Representation	Win/Mac		•		•
HELM Library Improvements	Win/Mac		•	•	•
Support for Ambiguous FASTA/HELM Monomers	Win/Mac		•		•
DNA/RNA Sequence Numbering Improvements	Win/Mac		•		•
Ring Fill Color Transfer to 3MF model **	Win/Mac				
3MF Export Settings for 3D Printing	Win/Mac			•	•
pic2mol Image to Structure Add-in (requires pic2mol subscription)	Win/Mac				
ChemOffice+ ***	Win/Mac				
Dedicated Tenant in multi-tenant Cloud Environment	Win/Mac				
ChemDraw JS for internal development	Win/Mac				
ChemOffice+ Features **	Platform	ChemDraw Prime	ChemDraw Professional	ChemOffice	ChemOffice+ Cloud Standard
Browse & Drill-down into ChemDraw Files (.cdx, .cdxml)	Win/Mac				
Browse & Drill-down ChemDraw Files embedded in MS Word	Win/Mac				•
Browse & Drill-down ChemDraw Files embedded in MS Powerpoint	Win/Mac				•
Browse ChemDraw For Excel Files	Win				•
Create Collection from .csv files	Win/Mac				
Create collection from SMILES text file	Win/Mac				•
Browse .mol & .sdf Files	Win/Mac				
View .sdf Files properties	Win/Mac				•
Copy Embedded Chemical Structures to the Clipboard	Win/Mac				
Create Collection of Chemical Structures	Win/Mac				
Adding Properties to Collections	Win/Mac				•
Editing Properties of Collections	Win/Mac				•
Saving Collection Layout as a Template	Win/Mac				•
Batch-Editing of Multiple Chemical Structures in Collections	Win/Mac				•
Structure-searching inside Cloud-hosted MS Office documents	Win/Mac				•
Searching across Signals Notebook (SNB) Experiments* ***	Win/Mac				
Create Collection of Reactions from SNB Experiments	Win/Mac				•
Export Collections to SD Files (v2000, v3000)	Win/Mac				

ChemOffice+ Features ** Continued	Platform	ChemDraw Prime	ChemDraw Professional	ChemOffice	ChemOffice+ Cloud Standard
reate Powerpoint Reaction Report Slide from SNB Experiments ****	Win/Mac				
Create Powerpoint Molecule Report Slide from Collection	Win/Mac				•
Recent Additions	Platform	ChemDraw Prime	ChemDraw Professional	ChemOffice	ChemOffice+ Cloud Standard
Magic Hotkeys Enhancements	Win/Mac	•	•	•	•
Shortcuts Enhancements	Win/Mac	•	•		•
Join function improvements	Win/Mac	•	•	•	•
Smart Copy/Paste (SMILES, InChl, HELM)	Win/Mac	•			
Aromatic Cycle Display Toggle and Preferences	Win/Mac	•	•		•
Stereochemistry handling improvements	Win/Mac	•			
Improved Polymer Brackets (Average MW)	Win/Mac	•	•	•	•
Atom/Bond Color Highlighting	Win/Mac		•		
Ring-Fill Coloring	Win/Mac		•	•	•
Search into CAS SciFinder-n	Win/Mac		•		
Search into Reaxys	Win/Mac		•		
Improved HELM Monomer Toolbar	Win/Mac				
HELM Monomer Versioning Support	Win/Mac		•		•
Copy as 3D-printable object (.3MF)**	Win/Mac				
Atom/Bond Color Highlight Transfer to 3MF	Win/Mac			•	
Google Patents/Scholar Add-in	Win/Mac				
PubChem GHS Safety Add-in	Win/Mac			•	•
ChemDraw Add-ins SDK	Win/Mac				•
ChemDraw Add-ins Dynamic Download	Win/Mac			•	•
Support for Add-ins Token-based Authentication	Win/Mac				
Shared HELM Libraries	Win/Mac			•	•
Includes	Platform	ChemDraw Prime	ChemDraw Professional	ChemOffice	ChemOffice+ Cloud Standard
Read and Save as .cdx / .cdxml Files	Win/Mac	•			
Read and Save as .rxn Files (v2000, v3000)	Win/Mac	•			•
Read and Save as .skc Files	Win/Mac	•			
Read and Save as .mol Files (v2000, v3000)	Win/Mac	•	•	•	•
Read and Save as .sdf Files (v2000, v3000)	Win/Mac	•	•	•	•
Read and Save as .rdf Files (v2000, v3000)	Win/Mac	•	•	•	•
Save ChemDraw Style Sheet	Win/Mac	•	•	•	•
Structure Clean-up	Win/Mac	•		•	
Reaction Clean-up	Win/Mac	•	•	•	•

Includes	Platform	ChemDraw Prime	ChemDraw Professional	ChemOffice	ChemOffice+ Cloud Standard
Magic Hotkeys	Win/Mac	•	•	•	•
Chemical Bonds Tools	Win/Mac	•		•	•
Text Tool	Win/Mac	•	•	•	•
3D Perspective Tool	Win/Mac	•		•	
Chemical Rings Tools	Win/Mac	•	•	•	•
Arrow Tool	Win/Mac	•	•	•	
Orbitals Tool	Win/Mac	•		•	
Brackets Tool	Win/Mac	•	•	•	
Pen Tools	Win/Mac	•		•	•
Shapes Tool	Win/Mac	•	•	•	•
Chemical Polymers Tools	Win/Mac	•		•	•
Mass Fragmentation Tools	Win/Mac	•		•	•
Thin Layer Chromatography Tool	Win/Mac	•	•	•	•
Gel Electrophoresis Tool	Win/Mac	•		•	
Insert OLE Object in ChemDraw	Win	•	•	•	•
Copy ChemDraw Structures as OLE Object	Win	•		•	
Show Stereochemistry	Win/Mac	•	•	•	•
Relative Stereochemistry (ISIS compatibility)	Win/Mac	•		•	•
Reaction Interpretation	Win/Mac	•	•	•	•
Reaction Mapping	Win/Mac	•		•	•
Calculate MW	Win/Mac	•	•	•	•
Calculate Exact Mass	Win/Mac	•		•	•
Calculate Chemical Formula	Win/Mac	•	•	•	•
Calculate Elemental Analysis	Win/Mac	•		•	
Calculate m/z	Win/Mac	•	•	•	•
Copy/Paste as CDXML	Win/Mac	•	•	•	•
Copy/Paste as SMILES	Win/Mac	•	•	•	•
Copy/Paste as SYBYL (SLN)	Win/Mac	•		•	
Copy/Paste as InChI	Win/Mac	•	•	•	•
Copy/Paste as Mol File / Mol3000	Win/Mac	•		•	
pKa / Log P / Log S	Win/Mac	•	•	•	
Atom List Generic Structures (Enumeration)	Win/Mac	•		•	•
tPSA	Win/Mac	•	•	•	•
Variable Attachment Generic Structures (Enumeration)	Win/Mac	•	•	•	•
Label Repeating Units Generic Structures (Enumeration)	Win/Mac	•	•	•	•
Polymer Repeating Units Generic Structures (Enumeration)	Win/Mac	•	•	•	
Chemical Structures Templates	Win/Mac	•	•	•	•
Laboratory Equipment Templates	Win/Mac	•	•	•	•

Includes	Platform	ChemDraw Prime	ChemDraw Professional	ChemOffice	ChemOffice+ Cloud Standard
Analyze/Check Structures	Win/Mac	•	•	•	•
Expand/Contract Labels	Win/Mac	•	•	•	•
Define/Use Nicknames	Win/Mac	•	•	•	•
Document Metadata/Tagging	Win/Mac	•	•	•	•
MulCple ChemDraw Items Folder	Win/Mac	•	•	•	•
Multicenter Attachments	Win/Mac	•	•	•	•
Save as JPEG image	Win/Mac	•	•	•	•
Save as PNG image	Win/Mac	•	•	•	•
Save as TIFF image	Win/Mac	•	•	•	•
Save as Scalable Vector Graphics (SVG)	Win/Mac	•	•	•	•
Save as Encapsulated Post Script (EPS)	Win/Mac	•	•	•	•
Name-to-Structure / Structure-to-Name	Win/Mac		•	•	•
Predict 1H NMR	Win/Mac		•	•	•
Predict 13C NMR	Win/Mac		•	•	•
Search SciFinder	Win/Mac		•	•	•
Search SciFinder-n	Win/Mac		•		
Search Reaxys	Win/Mac		•	•	•
Reaction Stoichiometry Grid	Win/Mac		•		•
R-Group Table Generic Structures (Enumeration)	Win/Mac		•	•	•
BioDraw Toolbar	Win/Mac		•		•
cLogP	Win/Mac		•	•	•
HELM Toolbar	Win/Mac		•		•
Copy/Paste as HELM	Win/Mac		•		•
Copy/Paste as FASTA Peptide	Win/Mac		•	•	•
Copy/Paste as FASTA DNA/RNA	Win/Mac		•	•	•
Support for HELM notation	Win/Mac		•		•
CAS RN to Structure from ChemACX.com	Win/Mac		•	•	•
Enhanced Stereochemistry Support	Win/Mac		•	•	•
ChemDraw for Excel	Win		•	•	•
ChemDraw for Excel	Win		•		•
CombiChem for Excel	Win		•	•	•
Name-to-Structure / Structure-to-Name for ChemDraw for Excel	Win		•	•	•
Chem3D Professional	Win		•	•	•
ChemFinder Standard	Win		•	•	•
ChemScript	Win		•	•	•
PubChem GHS Safety Add-in	Win/Mac				•
Google Scholar / Google Patents Add-in	Win/Mac			•	•
Save as 3D-printable object (.3MF)	Win/Mac			•	

Includes	Platform	ChemDraw Prime	ChemDraw Professional	ChemOffice	ChemOffice+ Cloud Standard
Copy as 3D-printable object (.3MF)	Win/Mac			•	•
Transfer Atom/Bond Color Highlights to 3D-printable object	Win/Mac			•	•
ChemACX Explorer	Win/Mac			•	
Custom ChemDraw Add-ins SDK	Win/Mac			•	•
Support for Token-based Authentication of Add-ins	Win/Mac				
Shared HELM Libraries	Win/Mac			•	•
Mnova ChemDraw Edition	Win/Mac			•	•
Chem3D Ultra	Win			•	•
Chem3D Interface to Conflex	Win			•	•
Chem3D Interface to Autodock	Win			•	•
Chem3D Interface to GAMESS 2020	Win			•	•
Chem3D Interface to Gaussian 16W	Win			•	•
Chem3D Interface to MOPAC 2016	Win			•	•
ChemFinder Ultra	Win			•	•
ChemFinder for Oracle	Win			•	•
Explorer Window View in ChemFinder Ultra	Win			•	•
BioViz in ChemFinder Ultra	Win			•	•
Compound Profiles in ChemDraw Finder Ultra	Win			•	•
Clustering in ChemFinder Ultra	Win			•	•
Combine ChemFinder Query Hit Lists	Win			•	
ChemFinder Exports to MS Word/Excel	Win			•	•
ClogP/CMR for ChemDraw for Excel/Chem3D	Win		•	•	•
Molecular Networks (pKa/Log P/Log S) for ChemDraw for Excel/Chem3D	Win		•	•	•
Molecular Topology for ChemDraw for Excel/Chem3D	Win		•	•	•
ChemProp Std Properties for ChemDraw for Excel/Chem3D	Win		•	•	•
ChemProp Pro Properties for ChemDraw for Excel/Chem3D	Win			•	•

^{*} Requires a minimum of 5 users









^{**} Pasting a .3MF object in MS Powerpoint requires a MS Office 365 subscription

^{***} ChemOffice+ is a Cloud-Native application and gets updated automatically on a quarterly basis

^{****} Requires a Signals Notebook Account